Citra-Lock™

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**Katheter  
Lock  
Lösung**

**Anti Microbiel  
Anti Koagulanz  
Anti Biofilm**

**30% & 46,7% TSC  
(Trisodium Citrat)**

**Dirinco®**  
Always innovating

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**Study Aim:** To determine whether locking with citrate 30% could increase overall catheter patency and reduce catheter related infections in comparison to heparin. The secondary aim was to determine whether a difference in major bleeding episodes could be demonstrated between citrate and heparin.  
**Study Design:** Multicentre (10 centres), double blind, randomised controlled, 291 catheters (143 heparin and 148 citrate 30%), 98 tunnelled catheters, 193 non-tunnelled catheters from catheter days 16.547.

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

### No. of vials of Urokinase used per month <sup>1</sup>

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#### REFERENCES:

1. Ash SR, Concentrated Sodium Citrate (23%) for Catheter Lock. Haemodialysis International 4:22-31, 2000
2. Weijmer JC, Randomized, Clinical Trial Comparison of Trisodium Citrate 30% and Heparin as Catheter-Locking Solution in Haemodialysis Patients. J Am Soc Nephrol. Sep;16(9):2769-77, 2005
3. Weijmer JC, Superior Antimicrobial Activity of Trisodium Citrate over Heparin for Catheter Locking. Nephrol Dial Transplant 17: 2189-2195, 2002
4. Shanks RM, Catheter Lock Solutions Influence Staphylococcal Biofilm Formation on Abiotic Surfaces. Nephrol Dial Transplant 21, 2247-2255, 2006
5. Nolan JP, Reducing Catheter Related Bacteraemia in Haemodialysis. Vascular Access Soc. 5th Int. Congress of Vascular Access Soc, Nice, 2007
6. Bosma JW, Siegert CE, Peerbooms PG, Weijmer MC, Reduction of biofilm formation with trisodium citrate in haemodialysis catheters: a randomized controlled trial. Nephrol Dial Transplant, 2009

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Dirinco AG - Bern - Switzerland